

MEDFIX[®]
INTERNATIONAL, LLC

Reducing cost without compromising quality.

- Medfix Founder

Medfix International, LLC has been delivering quality, low cost instrumentation around the world for many years. As a result of our success in fulfilling our mission of reducing cost without compromising quality, we have been able to partner with the biggest suppliers in the United States and Germany for the distribution of general surgical instrumentation. With over 14,000 items in stock and an expanding catalog with thousands of instruments, we can outfit anything from a single suture removal station to the largest multi-discipline operating room anywhere in the world, from simple forceps to procedure-specific instrumentation trays.

We can cross-reference almost any surgical instrument in the world by company and part number.

The vast majority of our orders are shipped within 3 business days. If you want to reduce costs, improve lead times, and remove guesswork from outfitting your surgical service, contact us and find out how easy it is to do business with a leader in the industry.

QUALITY ASSURANCE

Our instruments are manufactured from the highest quality German surgical stainless steel (ASTM F899), are treated to be corrosion resistant, and meet or exceed the most stringent of FDA and CE certifications.

Our instruments have been tested with the Thermo Scientific™ Niton XL3t XRF analyzer and show no significant variation from industry standard for type 420 stainless steel.

GUARANTEE

We offer the broadest, highest quality line of German, hand-held, reusable surgical instruments. Because of that, our instruments have a lifetime warranty against defects in material or workmanship when used for their intended purpose and when cared for according to industry standards. If you have any questions about use or care please do not hesitate to call us.

RETURNS

We strive to make our return process simple. Call (866-321-6331 or 520-398-5467) or email us (sales@medfix.com) to get a Return Material Authorization (RMA) number and shipping instructions. Be ready with the reason for the return, with your contact information, and invoice number (if available). All returns must have an RMA number or the claim cannot be processed. Replacements are shipped as soon as possible (depending upon availability), usually within 3 days of receipt of the returned item. Please package your returned items securely to avoid shipping damage.

We cannot accept instruments for replacement which have been altered or mishandled by the customer.

OUR STORY

Success is a journey, not a destination.

Since its inception, the driving force behind Medfix International, LLC has been to be part of the solution, and not perpetuate the status quo.

We came to exist almost out of necessity and as an answer to the overwhelming demand for lower cost alternatives to the traditional sources for surgical equipment, devices, and supplies. As formidably large corporations continued to grow exponentially, so did their prices. While many positive strides were made in the delivery of cutting edge surgical technology, the companies which spearheaded these advancements also lost sight of the responsibility to ensure their customers (surgeons, hospitals, surgery centers) were able to keep their doors open. With pricing that would eventually become cost-prohibitive and, in some cases, force their customers completely out of business, these corporations allowed profits to become paramount and fiscal responsibility slowly faded into the background.

Our mission was clear—we would seek, find, and provide to our customers alternatives for the equipment, devices, and supplies they needed. We knew quality must not be compromised but could not lose sight of the importance of affordable pricing.

Our philosophy today remains unchanged from the early years. We exist to serve our global community. From a makeshift operating room in a small village in Nepal, to a 10 bed hospital in rural Colombia, to the most prestigious centers of medical excellence in the world, Medfix International, LLC is committed to deliver surgical instruments, procedure-specific instrumentation, and implants in the fields of our expertise, namely, orthopedic, spine, and neurosurgery.

Our passion lies in assisting our customers develop, implement, and improve their surgical service, always making sure the patient comes first, and profits come second. We have seen how the delivery of our solutions has enabled surgeons to perform surgical procedures their hospitals and patients could not have otherwise afforded. This philosophy has propelled us forward, earned the trust of thousands of customers worldwide, and made it possible to improve patient's lives.

To our loyal customers, we thank you for your continued trust and patronage, and we invite those who seek an alternative who will work with you and help you succeed, come join us on our journey to make a difference and do right by our respective communities.

ELEVATORS

Woodson Elevator and Spatula	page 108
Cushing Elevator	page 108
Cobb Spinal Elevator	page 109
Cobb Elevator	page 109
Langenbeck Elevator	page 109
Penfield Dissector	page 110
Cobb Elevator	page 110
Freer Elevator	page 110



MF225-2832

WOODSON
ELEVATOR
AND SPATULA

Cat. No.	Description
MF225-2832	255 mm / 10 inch / 3 mm double ended / blunt



MF225-2872

ADSON
ELEVATOR

Cat No.	Description
MF225-2872	170 mm / 6¾ inch / sharp straight edge / straight / 8 mm blade
MF225-2887	170 mm / 6¾ inch / semi-sharp curved edge / straight / 8 mm blade
MF225-2882	170 mm / 6¾ inch / semi-sharp curved edge / curved / 6 mm blade
MF225-2877	170 mm / 6¾ inch / blunt curved edge / slightly curved / 6 mm blade



MF225-2987

CUSHING
ELEVATOR

Cat No.	Description
MF225-2997	190 mm / 7½ inch / sharp / straight edge / 14 mm
MF225-2987	190 mm / 7½ inch / sharp / curved edge / 14.5 mm
MF225-2992	190 mm / 7½ inch / sharp / serrated / curved edge / 14.5 mm



MF225-3147

**COBB
SPINAL
ELEVATOR**

Cat No.	Description
MF225-3142	240 mm / 9½ inch / hollow handle / 10 mm / ⅜ inch wide
MF225-3147	240 mm / 9½ inch / hollow handle / 13 mm / ½ inch wide
MF225-3152	240 mm / 9½ inch / hollow handle / 19 mm / ¾ inch wide



MF225-3112

**COBB
ELEVATOR**

Cat No.	Description
MF225-3102	240 mm / 9½ inch / solid stainless hexagonal handle / 10 mm / ⅜ inch wide / sharp
MF225-3107	240 mm / 9½ inch / solid stainless hexagonal handle / 13 mm / ½ inch wide / sharp
MF225-3112	240 mm / 9½ inch / solid stainless hexagonal handle / 19 mm / ¾ inch wide / sharp
MF225-3117	240 mm / 9½ inch / solid stainless hexagonal handle / 25 mm / 1 inch wide / sharp
MF225-3122	240 mm / 9½ inch / solid stainless hexagonal handle / 32 mm / 1¼ inch wide / sharp



MF210-9507

**LANGENBECK
ELEVATOR**

Cat. No.	Description
MF210-9502	205 mm / 8 inch / 6 mm / blunt
MF210-9507	205 mm / 8 inch / 17 mm / sharp

Elevators



MF225-4382

PENFIELD
DISSECTOR

Cat No.	Description
MF225-4382	180 mm / 7 inch / #1 / double ended
MF225-4387	185 mm / 7¼ inch / #2 / double ended
MF225-4392	180 mm / 7 inch / #3 / double ended
MF225-4397	220 mm / 8¾ inch / #4 / single ended
MF225-4402	290 mm / 11½ inch / #5 / double ended



MF225-3412

COBB
ELEVATOR

Cat No.	Description
MF225-3402	280 mm / 11 inch / knurled stainless steel handle / sharp / 10 mm / ¾ inch wide
MF225-3407	280 mm / 11 inch / knurled stainless steel handle / sharp / 13 mm / ½ inch wide
MF225-3412	280 mm / 11 inch / knurled stainless steel handle / sharp / 19 mm / ¾ inch wide
MF225-3417	280 mm / 11 inch / knurled stainless steel handle / sharp / 25 mm / 1 inch wide
MF225-3422	280 mm / 11 inch / knurled stainless steel handle / sharp / 32 mm / 1¼ inch wide



MF303-3402

FREER
ELEVATOR

Cat. No.	Description
MF303-3402	190 mm / 7½ inch / side wings prevent rolling / sharp and blunt 5 mm ends

NOTES

This image shows a full page of blank, lined paper. It features approximately 20 evenly spaced horizontal grey lines across its entire width, typical of notebook or legal stationery. The background is a solid off-white color, and there are no margins, text, or other markings present.

NOTES

[illegible]