

**MEDFIX<sup>®</sup>**  
INTERNATIONAL, LLC

Reducing cost without  
compromising quality.

*- Medfix Founder*

Medfix International, LLC has been delivering quality, low cost instrumentation around the world for many years. As a result of our success in fulfilling our mission of reducing cost without compromising quality, we have been able to partner with the biggest suppliers in the United States and Germany for the distribution of general surgical instrumentation. With over 14,000 items in stock and an expanding catalog with thousands of instruments, we can outfit anything from a single suture removal station to the largest multi-discipline operating room anywhere in the world, from simple forceps to procedure-specific instrumentation trays.

We can cross-reference almost any surgical instrument in the world by company and part number.

The vast majority of our orders are shipped within 3 business days. If you want to reduce costs, improve lead times, and remove guesswork from outfitting your surgical service, contact us and find out how easy it is to do business with a leader in the industry.

#### QUALITY ASSURANCE

Our instruments are manufactured from the highest quality German surgical stainless steel (ASTM F899), are treated to be corrosion resistant, and meet or exceed the most stringent of FDA and CE certifications.

Our instruments have been tested with the Thermo Scientific™ Niton XL3t XRF analyzer and show no significant variation from industry standard for type 420 stainless steel.

#### GUARANTEE

**We offer the broadest, highest quality line of German, hand-held, reusable surgical instruments. Because of that, our instruments have a lifetime warranty against defects in material or workmanship when used for their intended purpose and when cared for according to industry standards. If you have any questions about use or care please do not hesitate to call us.**

#### RETURNS

We strive to make our return process simple. Call (866-321-6331 or 520-398-5467) or email us (sales@medfix.com) to get a Return Material Authorization (RMA) number and shipping instructions. Be ready with the reason for the return, with your contact information, and invoice number (if available). All returns must have an RMA number or the claim cannot be processed. Replacements are shipped as soon as possible (depending upon availability), usually within 3 days of receipt of the returned item. Please package your returned items securely to avoid shipping damage.

We cannot accept instruments for replacement which have been altered or mishandled by the customer.

## OUR STORY

Success is a journey, not a destination.

Since its inception, the driving force behind Medfix International, LLC has been to be part of the solution, and not perpetuate the status quo.

We came to exist almost out of necessity and as an answer to the overwhelming demand for lower cost alternatives to the traditional sources for surgical equipment, devices, and supplies. As formidably large corporations continued to grow exponentially, so did their prices. While many positive strides were made in the delivery of cutting edge surgical technology, the companies which spearheaded these advancements also lost sight of the responsibility to ensure their customers (surgeons, hospitals, surgery centers) were able to keep their doors open. With pricing that would eventually become cost-prohibitive and, in some cases, force their customers completely out of business, these corporations allowed profits to become paramount and fiscal responsibility slowly faded into the background.

Our mission was clear—we would seek, find, and provide to our customers alternatives for the equipment, devices, and supplies they needed. We knew quality must not be compromised but could not lose sight of the importance of affordable pricing.


Our philosophy today remains unchanged from the early years. We exist to serve our global community. From a makeshift operating room in a small village in Nepal, to a 10 bed hospital in rural Colombia, to the most prestigious centers of medical excellence in the world, Medfix International, LLC is committed to deliver surgical instruments, procedure-specific instrumentation, and implants in the fields of our expertise, namely, orthopedic, spine, and neurosurgery.

Our passion lies in assisting our customers develop, implement, and improve their surgical service, always making sure the patient comes first, and profits come second. We have seen how the delivery of our solutions has enabled surgeons to perform surgical procedures their hospitals and patients could not have otherwise afforded. This philosophy has propelled us forward, earned the trust of thousands of customers worldwide, and made it possible to improve patient's lives.

To our loyal customers, we thank you for your continued trust and patronage, and we invite those who seek an alternative who will work with you and help you succeed, come join us on our journey to make a difference and do right by our respective communities.

## **KERRISONS/RONGEURS**

Kerrison Swing Rongeur	page 18
Micro Kerrison 40° up-cutting	page 19
Micro Kerrison 90° up-cutting	page 20
Micro Kerrison 40° down-cutting	page 21
Micro Kerrison 90° down-cutting	page 22
Spurling Kerrison	page 23
Bayonet Kerrison	page 24
Foraminotomy Kerrison	page 25
Bayonet Foraminotomy Kerrison	page 26
Kerrison Side-Bayonet	page 27
Cushing IVD Rongeur	page 28
Bayonet Cushing IVD Rongeur	page 28
Love Gruenwald IVD Rongeur	page 29
Bayonet Love Gruenwald Rongeur	page 29
Spurling IVD Rongeur	page 30
Bayonet Spurling IVD Rongeur	page 30
Caspar Rongeur	page 31
Caspar Fenestrated Serrated Rongeur	page 32
Ferris-Smith Spurling Rongeur	page 33
Decker IVD Rongeur	page 34
Peapod IVD Rongeur	page 35



STANDARD FINISH  
SWING

SURFACE / Stainless Steel  
HOW TO ORDER / MF249-5101



BLACK FINISH  
SWING

SURFACE / Black / ceramic coated  
HOW TO ORDER / MF249-5101CC



STANDARD THIN FOOTPLATE

SURFACE / Stainless Steel  
HOW TO ORDER / MF249-5101T



BLACK THIN FOOTPLATE

SURFACE / Black / ceramic coated  
HOW TO ORDER / MF245-7202TCC



MF249-7201

MEDFIX®  
**KERRISON**  
SWING RONGEUR

	180 mm / 7 inch	200 mm / 8 inch	Description
 1 mm	MF249-7101	MF249-7201	40° up-cutting
 2 mm	MF249-7102	MF249-7202	40° up-cutting
 3 mm	MF249-7103	MF249-7203	40° up-cutting
 4 mm	MF249-7104	MF249-7204	40° up-cutting
 5 mm	MF249-7105	MF249-7205	40° up-cutting
 6 mm	MF249-7106	MF249-7206	40° up-cutting
 1 mm	MF249-5101	MF249-5201	90° up-cutting
 2 mm	MF249-5102	MF249-5202	90° up-cutting
 3 mm	MF249-5103	MF249-5203	90° up-cutting
 4 mm	MF249-5104	MF249-5204	90° up-cutting
 5 mm	MF249-5105	MF249-5205	90° up-cutting
 6 mm	MF249-5106	MF249-5206	90° up-cutting



STANDARD FINISH SWING  
SURFACE / Stainless Steel  
HOW TO ORDER / MF249-5101



BLACK FINISH SWING  
SURFACE / Black / ceramic coated  
HOW TO ORDER / MF249-5101CC



STANDARD THIN FOOTPLATE  
SURFACE / Stainless Steel  
HOW TO ORDER / MF249-5101T



BLACK THIN FOOTPLATE  
SURFACE / Black / ceramic coated  
HOW TO ORDER / MF245-7202TCC



MEDFIX®  
MICRO  
KERRISON  
with ejector

small handle

180 mm / 7 inch

200 mm / 8 inch

230 mm / 9 inch

Description

-  1 mm
-  2 mm
-  3 mm
-  4 mm
-  5 mm
-  6 mm

MF235-7201

MF245-7201

MF239-9001

40° up-cutting

MF235-7202

MF245-7202

MF239-9002

40° up-cutting

MF235-7203

MF245-7203

MF239-9003

40° up-cutting

MF235-7204

MF245-7204

MF239-9004

40° up-cutting

MF235-7205

MF245-7205

MF239-9005

40° up-cutting

MF235-7206

MF245-7206

MF239-9006

40° up-cutting

STANDARD FINISH



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B

STANDARD THIN FOOTPLATE



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202T

BLACK THIN FOOTPLATE



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202TB



Up-cutting / 90°

MF290-7001

**MEDFIX®  
MICRO  
KERRISON  
with ejector**

*small handle*

180 mm / 7 inch

200 mm / 8 inch

230 mm / 9 inch

Description

-  1 mm
-  2 mm
-  3 mm
-  4 mm
-  5 mm
-  6 mm

MF290-7001  
MF290-7002  
MF290-7003  
MF290-7004  
MF290-7005  
MF290-7006

MF290-8001  
MF290-8002  
MF290-8003  
MF290-8004  
MF290-8005  
MF290-8006

MF290-9001  
MF290-9002  
MF290-9003  
MF290-9004  
MF290-9005  
MF290-9006

90° up-cutting  
90° up-cutting  
90° up-cutting  
90° up-cutting  
90° up-cutting  
90° up-cutting

STANDARD FINISH



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B

STANDARD THIN FOOTPLATE



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202T

BLACK THIN FOOTPLATE



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202TB



Down-cutting / 40°

MF241-7403

**MEDFIX®  
MICRO  
KERRISON  
with ejector**

small handle	180 mm / 7 inch	200 mm / 8 inch	230 mm / 9 inch	Description
 1 mm	MF241-7301	MF241-7401	MF241-7501	40° down-cutting
 2 mm	MF241-7302	MF241-7402	MF241-7502	40° down-cutting
 3 mm	MF241-7303	MF241-7403	MF241-7503	40° down-cutting
 4 mm	MF241-7304	MF241-7404	MF241-9004	40° down-cutting
 5 mm	MF241-7305	MF241-7405	MF241-7505	40° down-cutting
 6 mm	MF241-7306	MF241-7406	MF241-7506	40° down-cutting

For additional options please see the index under Micro Kerrison.

STANDARD FINISH 

SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH 

SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B

STANDARD THIN FOOTPLATE 

SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202T

BLACK THIN FOOTPLATE 

SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202TB



Down-cutting / 90°

MF238-7301

MEDFIX®  
MICRO  
**KERRISON**  
with ejector

small handle

180 mm / 7 inch

200 mm / 8 inch

230 mm / 9 inch

Description

-  1 mm
-  2 mm
-  3 mm
-  4 mm
-  5 mm
-  6 mm

MF238-7301  
MF238-7302  
MF238-7303  
MF238-7304  
MF238-7305  
MF238-7306

MF238-7401  
MF238-7402  
MF238-7403  
MF238-7404  
MF238-7405  
MF238-7406

MF238-7501  
MF238-7502  
MF238-7503  
MF238-7504  
MF238-7505  
MF238-7506

90° down-cutting  
90° down-cutting  
90° down-cutting  
90° down-cutting  
90° down-cutting  
90° down-cutting

STANDARD FINISH



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B

STANDARD THIN FOOTPLATE



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202T

BLACK THIN FOOTPLATE



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202TB




MF231-3907

black ceramic coated

**MEDFIX®  
SPURLING  
KERRISON**  
thin footplate

	155 mm / 6 inch	180 mm / 7 inch	200 mm / 8 inch	200 mm / 8 inch	Description
2 mm	MF231-3702	MF231-3867	MF231-4027	MF232-3027B	90° up-cutting / 15 mm opening
3 mm	MF231-3707	MF231-3872	MF231-4032	MF232-3032B	90° up-cutting / 15 mm opening
4 mm	MF231-3712	MF231-3877	MF231-4037	MF232-3037B	90° up-cutting / 15 mm opening
5 mm	MF231-3717	MF231-3882	MF231-4042	MF232-3042B	90° up-cutting / 15 mm opening
6 mm	MF231-3722	MF231-3887	MF231-4047	MF232-3047B	90° up-cutting / 15 mm opening
2 mm	MF231-3742	MF231-3907	MF231-4067		90° up-cutting / 15 mm opening
3 mm	MF231-3747	MF231-3912	MF231-4072		90° up-cutting / 15 mm opening
4 mm	MF231-3752	MF231-3917	MF231-4077		90° up-cutting / 15 mm opening
5 mm	MF231-3757	MF231-3922	MF231-4082		90° up-cutting / 15 mm opening
6 mm	MF231-3762	MF231-3927	MF231-4087		90° up-cutting / 15 mm opening
2 mm	MF231-3822	MF231-3987	MF231-4142		90° down-cutting / 15 mm opening
3 mm	MF231-3827	MF231-3992	MF231-4147		90° down-cutting / 15 mm opening
4 mm	MF231-3832	MF231-3997	MF231-4152		90° down-cutting / 15 mm opening
5 mm	MF231-3837	MF231-4002	MF231-4157		90° down-cutting / 15 mm opening
6 mm	MF231-3842	MF231-4007	MF231-4162		90° down-cutting / 15 mm opening
	230 mm / 9 inch 40° up-cutting	230 mm / 9 inch 90° up-cutting	280 mm / 11 inch 40° up-cutting	Description	
2 mm	MF231-4182	MF231-4222	MF231-4507	large handle / 15 mm opening	
3 mm	MF231-4187	MF231-4227	MF231-4512	large handle / 15 mm opening	
4 mm	MF231-4192	MF231-4232	MF231-4517	large handle / 15 mm opening	
5 mm	MF231-4197	MF231-4237	MF231-4522	large handle / 15 mm opening	
6 mm	MF231-4202	MF231-4242	MF231-4527	large handle / 15 mm opening	



MF367-3602B

MEDFIX®  
BAYONET  
KERRISON  
with ejector

	110 mm / 4 1/8 inch	150 mm / 6 inch	180 mm / 7 inch	Description
 1 mm	MF267-1601	MF267-2601	MF267-3601	90° up-cutting
 2 mm	MF267-1602	MF267-2602	MF267-3602	90° up-cutting
 3 mm	MF267-1603	MF267-2603	MF267-3603	90° up-cutting
 4 mm	MF267-1604	MF267-2604	MF267-3604	90° up-cutting
 5 mm	MF267-1605	MF267-2605	MF267-3605	90° up-cutting
 6 mm	MF267-1606	MF267-2606	MF267-3606	90° up-cutting
 1 mm	MF367-1601	MF367-2601	MF367-3601	40° up-cutting
 2 mm	MF367-1602	MF367-2602	MF367-3602	40° up-cutting
 3 mm	MF367-1603	MF367-2603	MF367-3603	40° up-cutting
 4 mm	MF367-1604	MF367-2604	MF367-3604	40° up-cutting
 5 mm	MF367-1605	MF367-2605	MF367-3605	40° up-cutting
 6 mm	MF367-1606	MF367-2606	MF367-3606	40° up-cutting

STANDARD FINISH



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B



MEDFIX®  
FORAMINOTOMY  
**KERRISON**  
with ejector

	180 mm / 7 inch	200 mm / 8 inch	230 mm / 9 inch	Description
 1 mm	MF676-7501	MF676-8501	MF676-9501	130° curved up / 40° up-cutting
 2 mm	MF676-7502	MF676-8502	MF676-9502	130° curved up / 40° up-cutting
 3 mm	MF676-7503	MF676-8503	MF676-9503	130° curved up / 40° up-cutting
 4 mm	MF676-7504	MF676-8504	MF676-9504	130° curved up / 40° up-cutting
 5 mm	MF676-7505	MF676-8505	MF676-9505	130° curved up / 40° up-cutting
 6 mm	MF676-7506	MF676-8506	MF676-9506	130° curved up / 40° up-cutting
 1 mm	MF776-7501	MF776-8501	MF776-9501	130° curved down / 40° down-cutting
 2 mm	MF776-7502	MF776-8502	MF776-9502	130° curved down / 40° down-cutting
 3 mm	MF776-7503	MF776-8503	MF776-9503	130° curved down / 40° down-cutting
 4 mm	MF776-7504	MF776-8504	MF776-9504	130° curved down / 40° down-cutting
 5 mm	MF776-7505	MF776-8505	MF776-9505	130° curved down / 40° down-cutting
 6 mm	MF776-7506	MF776-8506	MF776-9506	130° curved down / 40° down-cutting

STANDARD FINISH




SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B



MF676-4001B

MEDFIX®  
BAYONET  
FORAMINOTOMY  
**KERRISON**  
with ejector

	110 mm / 4½ inch	150 mm / 6 inch	180 mm / 7 inch	Description
 1 mm	MF676-4001	MF676-6001	MF676-7001	130° curved up / 40° up-cutting
 2 mm	MF676-4002	MF676-6002	MF676-7002	130° curved up / 40° up-cutting
 3 mm	MF676-4003	MF676-6003	MF676-7003	130° curved up / 40° up-cutting
 4 mm	MF676-4004	MF676-6004	MF676-7004	130° curved up / 40° up-cutting
 5 mm	MF676-4005	MF676-6005	MF676-7005	130° curved up / 40° up-cutting
 6 mm	MF676-4006	MF676-6006	MF676-7006	130° curved up / 40° up-cutting
 1 mm	MF776-4001	MF776-6001	MF776-7001	130° curved down / 40° down-cutting
 2 mm	MF776-4002	MF776-6002	MF776-7002	130° curved down / 40° down-cutting
 3 mm	MF776-4003	MF776-6003	MF776-7003	130° curved down / 40° down-cutting
 4 mm	MF776-4004	MF776-6004	MF776-7004	130° curved down / 40° down-cutting
 5 mm	MF776-4005	MF776-6005	MF776-7005	130° curved down / 40° down-cutting
 6 mm	MF776-4006	MF776-6006	MF776-7006	130° curved down / 40° down-cutting

STANDARD FINISH 

SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH 

SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B




MEDFIX®  
**KERRISON**  
SIDE-BAYONET  
with ejector

	110 mm / 4 1/8 inch	150 mm / 6 inch	180 mm / 7 inch	Description
1 mm	MF476-2101	MF476-3101	MF476-4101	left / 40° up-cutting
2 mm	MF476-2102	MF476-3102	MF476-4102	left / 40° up-cutting
3 mm	MF476-2103	MF476-3103	MF476-4103	left / 40° up-cutting
4 mm	MF476-2104	MF476-3104	MF476-4104	left / 40° up-cutting
5 mm	MF476-2105	MF476-3105	MF476-4105	left / 40° up-cutting
6 mm	MF476-2106	MF476-3106	MF476-4106	left / 40° up-cutting
1 mm	MF276-2101	MF276-3101	MF276-4101	left / 90° up-cutting
2 mm	MF276-2102	MF276-3102	MF276-4102	left / 90° up-cutting
3 mm	MF276-2103	MF276-3103	MF276-4103	left / 90° up-cutting
4 mm	MF276-2104	MF276-3104	MF276-4104	left / 90° up-cutting
5 mm	MF276-2105	MF276-3105	MF276-4105	left / 90° up-cutting
6 mm	MF276-2106	MF276-3106	MF276-4106	left / 90° up-cutting
	110 mm / 4 1/8 inch	150 mm / 6 inch	180 mm / 7 inch	
1 mm	MF576-2101	MF576-3101	MF576-4101	right / 40° up-cutting
2 mm	MF576-2102	MF576-3102	MF576-4102	right / 40° up-cutting
3 mm	MF576-2103	MF576-3103	MF576-4103	right / 40° up-cutting
4 mm	MF576-2104	MF576-3104	MF576-4104	right / 40° up-cutting
5 mm	MF576-2105	MF576-3105	MF576-4105	right / 40° up-cutting
6 mm	MF576-2106	MF576-3106	MF576-4106	right / 40° up-cutting
1 mm	MF376-2101	MF376-3101	MF376-4101	right / 90° up-cutting
2 mm	MF376-2102	MF376-3102	MF376-4102	right / 90° up-cutting
3 mm	MF376-2103	MF376-3103	MF376-4103	right / 90° up-cutting
4 mm	MF376-2104	MF376-3104	MF376-4104	right / 90° up-cutting
5 mm	MF376-2105	MF376-3105	MF376-4105	right / 90° up-cutting
6 mm	MF376-2106	MF376-3106	MF376-4106	right / 90° up-cutting



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202


SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B



MF232-2347

**MEDFIX®**  
**CUSHING**  
IVD **RONGEUR**

Straight	Angled Up	Angled Down	Description
MF221-2001	MF221-2000	MF221-1999	180 mm / 7 inch / 1.5 x 10 mm bite
MF221-1998	MF221-1997	MF221-1996	230 mm / 9 inch / 1.5 x 10 mm bite
MF232-2302	MF232-2307	MF232-2312	130 mm / 5 inch / 2 x 10 mm bite
MF232-2322	MF232-2327	MF232-2332	150 mm / 6 inch / 2 x 10 mm bite
MF232-2347	MF232-2352	MF232-2357	180 mm / 7 inch / 2 x 10 mm bite
MF232-2367	MF232-2372	MF232-2377	230 mm / 9 inch / 2 x 10 mm bite
MF232-2387	MF232-2392	MF232-2394	255 mm / 10 inch / 2 x 10 mm bite
MF232-2412	MF232-2417	MF232-2422	280 mm / 11 inch / 2 x 10 mm bite
MF232-2427	MF232-2428	MF232-2329	305 mm / 12 inch / 2 x 10 mm bite



MF137-7202B

**MEDFIX®**  
**BAYONET**  
**CUSHING**  
IVD **RONGEUR**

Cat. No.	Description
MF137-7001	180 mm / 7 inch / 1.5 x 10 mm bite / straight
MF137-7101	180 mm / 7 inch / 1.5 x 10 mm bite / angled up
MF137-7201	180 mm / 7 inch / 1.5 x 10 mm bite / angled down
MF137-7002	180 mm / 7 inch / 2 x 10 mm bite / straight
MF137-7102	180 mm / 7 inch / 2 x 10 mm bite / angled up
MF137-7202	180 mm / 7 inch / 2 x 10 mm bite / angled down

STANDARD FINISH




SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH




SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B



MEDFIX®  
LOVE  
GRUENWALD  
IVD **RONGEUR**

Straight	Angled Up	Angled Down	Description
MF232-2432	MF232-2437	MF232-2442	130 mm / 5 inch / 3 x 10 mm bite
MF232-2452	MF232-2457	MF232-2462	150 mm / 6 inch / 3 x 10 mm bite
MF232-2477	MF232-2482	MF232-2487	180 mm / 7 inch / 3 x 10 mm bite
MF232-2497	MF232-2502	MF232-2507	230 mm / 9 inch / 3 x 10 mm bite
MF232-2517	MF232-2522	MF232-2527	255 mm / 10 inch / 3 x 10 mm bite
MF232-2537	MF232-2542	MF232-2572	280 mm / 11 inch / 3 x 10 mm bite
MF232-2552	MF232-2553	MF232-2554	305 mm / 12 inch / 3 x 10 mm bite



MEDFIX®  
BAYONET  
LOVE  
GRUENWALD  
**RONGEUR**


Straight	Angled Up	Angled Down	Description
MF138-7001	MF138-7002	MF138-7003	180 mm / 7 inch / 3 x 10 mm bite

STANDARD FINISH   
SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH   
SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B

STANDARD FINISH  
DETACHABLE   
SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202DT


BLACK FINISH  
DETACHABLE   
SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202DTB



MF232-2607

MEDFIX®  
SPURLING  
IVD RONGEUR

	Straight	Angled Up	Angled Down	Description
	MF232-2557	MF232-2562	MF232-2567	130 mm / 5 inch / 4 x 10 mm bite
	MF232-2577	MF232-2582	MF232-2587	150 mm / 6 inch / 4 x 10 mm bite
	MF232-2602	MF232-2607	MF232-2612	180 mm / 7 inch / 4 x 10 mm bite
	MF232-2622	MF232-2627	MF232-2632	230 mm / 9 inch / 4 x 10 mm bite
	MF232-2642	MF232-2647	MF232-2652	255 mm / 10 inch / 4 x 10 mm bite
	MF232-2662	MF232-2667	MF232-2668	280 mm / 11 inch / 4 x 10 mm bite
	MF232-2674	MF232-2676	MF232-2678	305 mm / 12 inch / 4 x 10 mm bite



MF139-7101B

MEDFIX®  
BAYONET  
SPURLING  
IVD  
RONGEUR

Cat. No.	Description
MF139-7001	180 mm / 7 inch / 4 x 10 mm bite / straight
MF139-7002	180 mm / 7 inch / 4 x 10 mm bite / angled up
MF139-7003	180 mm / 7 inch / 4 x 10 mm bite / angled down

STANDARD FINISH



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B



MF331-4002

**MEDFIX®  
CASPAR  
RONGEUR**

Straight	Angled Up	Angled Down	Description
MF331-4002	MF331-5002	MF331-6002	2 x 15 mm bite / 120 mm / 4¾ inch
MF332-4002	MF332-5002	MF332-6002	3 x 15 mm bite / 120 mm / 4¾ inch
MF333-4002	MF333-5002	MF333-6002	4 x 15 mm bite / 120 mm / 4¾ inch
MF334-4002	MF334-5002	MF334-6002	5 x 15 mm bite / 120 mm / 4¾ inch
MF336-4002	MF336-5002	MF336-6002	2 x 15 mm bite / 155 mm / 6 inch
MF336-4102	MF336-5102	MF336-6102	3 x 15 mm bite / 155 mm / 6 inch
MF336-4202	MF336-5202	MF336-6202	4 x 15 mm bite / 155 mm / 6 inch
MF336-4302	MF336-5302	MF336-6302	5 x 15 mm bite / 155 mm / 6 inch
MF337-4002	MF337-5002	MF337-6002	2 x 15 mm bite / 180 mm / 7 inch
MF337-4502	MF337-5502	MF337-6502	3 x 15 mm bite / 180 mm / 7 inch
MF337-4702	MF337-5702	MF337-6702	4 x 15 mm bite / 180 mm / 7 inch
MF337-4802	MF337-5802	MF337-6805	5 x 15 mm bite / 180 mm / 7 inch

STANDARD FINISH



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B




**MEDFIX®**  
**CASPAR**  
FENESTRATED  
**RONGEUR**  
Serrated

MF618-5801

Serrated	Straight	Angled Up	Angled Down	Description
	MF606-5201	MF606-5301	MF606-5401	2 x 12 mm bite / 140 mm / 5½ inch
	MF607-5201	MF607-5301	MF607-5401	3 x 12 mm bite / 140 mm / 5½ inch
	MF608-5201	MF608-5301	MF608-5401	4 x 14 mm bite / 140 mm / 5½ inch
	MF609-5201	MF609-5301	MF609-5401	5 x 14 mm bite / 140 mm / 5½ inch
	MF610-5201	MF610-5301	MF610-5401	6 x 16 mm bite / 140 mm / 5½ inch
	MF612-5201	MF613-5301	MF614-5401	2 x 12 mm bite / 160 mm / 6¼ inch
	MF615-5501	MF616-5601	MF617-5701	3 x 12 mm bite / 160 mm / 6¼ inch
	MF618-5801	MF618-5901	MF619-4001	4 x 14 mm bite / 160 mm / 6¼ inch
	MF619-5001	MF619-5002	MF619-5003	5 x 14 mm bite / 160 mm / 6¼ inch
	MF619-6001	MF619-6002	MF619-6003	6 x 16 mm bite / 160 mm / 6¼ inch
	MF620-7001	MF620-7002	MF620-7003	2 x 12 mm bite / 185 mm / 7¼ inch
	MF623-7001	MF623-7002	MF623-7003	3 x 12 mm bite / 185 mm / 7¼ inch
	MF624-7001	MF624-7002	MF624-7003	4 x 14 mm bite / 185 mm / 7¼ inch
	MF625-7001	MF625-7002	MF625-7003	5 x 14 mm bite / 185 mm / 7¼ inch
	MF626-7001	MF626-7002	MF626-7003	6 x 16 mm bite / 185 mm / 7¼ inch





MEDFIX®  
FERRIS-SMITH  
SPURLING  
**RONGEUR**

MF278-2201

Straight	Angled Up	Angled Down	Description
MF278-1801	MF278-1802	MF278-1803	180 mm / 7 inch / 2 x 10 mm bite
MF278-1901	MF278-1902	MF278-1903	180 mm / 7 inch / 3 x 10 mm bite
MF278-2001	MF278-2002	MF278-2003	180 mm / 7 inch / 4 x 10 mm bite
MF278-2101	MF278-2102	MF278-2103	180 mm / 7 inch / 5 x 10 mm bite
MF278-2101	MF278-2102	MF278-2103	180 mm / 7 inch / 5 x 10 mm bite
MF278-2201	MF278-2202	MF278-2203	180 mm / 7 inch / 6 x 10 mm bite

STANDARD FINISH




SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202

BLACK FINISH



SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B



MF271-5202

MEDFIX®  
DECKER  
IVD RONGEUR

Cat. No.	Description
MF271-5202	150 mm / 6 inch / 2 x 6 mm bite / straight
MF271-5207	150 mm / 6 inch / 2 x 6 mm bite / angled right
MF271-5212	150 mm / 6 inch / 2 x 6 mm bite / angled left
MF271-5217	150 mm / 6 inch / 2 x 6 mm bite / angled up
MF271-5222	150 mm / 6 inch / 2 x 6 mm bite / angled down
MF271-6097	150 mm / 6 inch / 1.5 x 6 mm bite / straight
MF272-5001*	140 mm / 5½ inch / 1.5 x 5.0 mm bite / straight
MF272-5002*	140 mm / 5½ inch / 1.5 x 5.0 mm bite / angled left
MF272-5003*	140 mm / 5½ inch / 1.5 x 5.0 mm bite / angled right
MF272-5004*	140 mm / 5½ inch / 1.5 x 5.0 mm bite / angled up
MF272-5005*	140 mm / 5½ inch / 1.5 x 5.0mm bite / angled down

\*THESE DECKER IVD RONGEURS DON'T COME WITH A BLACK FINISH OPTION.



STANDARD FINISH  
SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202



BLACK FINISH  
SURFACE / Black / titanium coated  
HOW TO ORDER / MF245-7202B



MF232-2807

MEDFIX®  
PEAPOD  
IVD **RONGEUR**

Cat. No.	Description
MF232-2791	180 mm / 7 inch / small / 2 mm bite
MF232-2800	130 mm / 5½ inch / extra small / 1.5 mm bite
MF232-2802	130 mm / 5½ inch / small / 2 mm bite
MF232-2803	180 mm / 7 inch / medium / 2.5 mm bite
MF232-2807	130 mm / 5½ inch / medium / 2.5 mm bite
MF232-2808	180 mm / 7 inch / large / 3.2 mm bite
MF232-2812	130 mm / 5½ inch / large / 3.2 mm bite

STANDARD FINISH



SURFACE / Stainless Steel  
HOW TO ORDER / MF245-7202





